

# Y and Tee Piece Extension Tubes in ASME BPE for Autoclave Chambers

## Product Code: YPE-BPE TPE-BPE





## **Description:**

Limited entry points on an autoclave chamber wall can present the problem of feeding both validation thermocouples and providing access to the chamber environment for a pressure transducer.

Extension tubes increase the flexibility of feeding thermocouples and pressure transducers through autoclave entry points.

ASME BPE Tee pieces are available in a range of standard sizes, see page 2

Y pieces and custom Tee pieces can also be manufactured on request.

Welding is carried out using orbital TIG welding.

Extension tubes can be supplied with full documentation, including certificate of conformance, material test certificates and 3.1B certificate.

#### Applications:

Pharmaceutical autoclave chambers.

Specification:			
MATERIAL	Tubing in 316L grade stainless steel to ASME BPE SF1 finish.		
FINISH	Highly polished with an inner surface of 0.5μm and outer surface finish 0.8μM.		
TUBE SIZE	<sup>3</sup> ⁄ <sub>4</sub> ", 1", 1½", 2", 2½" or 3".		
PROCESS CONNECTION	Tri-clover flanges to suit tube size, machined from round bar stock in 316L stainless steel to ASME BPE and prepared for orbital welding. Inner surface finished to Ra 0.5µM.		
FINAL FINISH	On request the completed tee / Y piece can be pickled, passivated and electropolished. This removes a layer of base material and improves the aesthetic appearance		

Datasheet: 035a




Fax: +44 (0) 1642 618307 Tel: +44 (0) 1642 602878 Email: sales@thermal-detection.com

Web: <u>www.thermal-detection.com</u>
© Thermal Detection Ltd 2016




# Y and Tee Piece Extension Tubes in **ASME BPE for Autoclave Chambers**

# **Equal and Short Outlet Tee Piece Standard Dimensions**Note: Custom sizes available on request



SIZE TUBE O.D.	A (mm)	B (mm)	C (mm)
1/2"	57.2	114.3	25.4
3/4"	60.3	120.6	28.6
1"	66.7	133.4	28.6
1½"	73	146	34.9
2"	85.7	171.5	41.3
2½"	92.1	184.2	47.6
3"	98.4	196.9	54



Datasheet: 035a

**Thermal Detection Limited** Unit 6, Primrose Hill Industrial Estate Orde Wingate Way Stockton-on-Tees, England, **TS19 0GA** 



Fax: +44 (0) 1642 618307 Tel: +44 (0) 1642 602878

Email: sales@thermal-detection.com Web: www.thermal-detection.com © Thermal Detection Ltd 2016



#### **Y and Tee Piece Extension Tubes in ASME BPE for Autoclave Chambers**

#### **Validation Accessories:**



Thermocouple Entry Gland TCEG



Pre-cut Validation Thermocouples



Validation Thermocouple Looms



Steriprobe Autoclave Load Probe



SL160 Angled Autoclave Probe



GL45-1 Bottle Cap Gland for Thermocouple



Clamps, Tri-clover Blanks, Gaskets



Type T thermocouple Wire



Ampoule Holders

#### Datasheet: 035a

**Thermal Detection Limited** Unit 6, Primrose Hill Industrial Estate Orde Wingate Way Stockton-on-Tees, England, **TS19 0GA** 



Fax: +44 (0) 1642 618307 Tel: +44 (0) 1642 602878 Email: sales@thermal-detection.com

Web: www.thermal-detection.com

© Thermal Detection Ltd 2016