

Standard Test Pack

Product Code: STP STPW STPC





Description:

The Standard Test Pack (STP) conforms to EN 285, EN 13060, CFPP 01-01 and HTM 2010 and is designed to enable the Bowie and Dick test to be carried out on large steam sterilisers.

The Standard Test Pack Wrap (STPW) is a useful accessory for the standard test pack. It helps the operator to ensure the sheets are correctly compressed and secure (EN 285 stipulates that the Standard Test Pack be compressed by hand, wrapped in a similar fabric and secured with tape that is no more than 25 mm wide). It is made with the same material as the standard test pack.

Ties and loops remove the need to use autoclave or masking tape to hold the pack together and allow an even compression of the pack. A handle allows the pack to be moved easily when transporting or placing in a sterilization chamber. Eyes are fitted at the top of the wrap to attach the free space temperature thermocouple aerial firmly in the correct position using the 50 mm height guide when performing small load thermometric tests.

The Standard Test Pack Carrier (STPC) provides protection for the Standard Test Pack and Standard Test Pack Wrap when being transported and whilst in storage.

Specification:	
MATERIAL	100% bleached cotton
SHEET SIZE	120cm x 90cm
PACK	38 sheets (35 required for test) including a thermocouple positioning sheet.
LAUNDRY	Service available as an extra for sheets to be prewashed, folded and ready for use.

Datasheet: 043

Thermal Detection Limited Unit 6, Primrose Hill Industrial Estate Orde Wingate Way Stockton-on-Tees, England, TS19 0GA



Fax: +44 (0) 1642 618307 Tel: +44 (0) 1642 602878 Email: sales@thermal-detection.com

Web: www.thermal-detection.com

© Thermal Detection Ltd 2016



Standard Test Pack

Validation Accessories:



Pre-cut Validation Thermocouples



Validation Thermocouple Looms



GL45 Bottle Cap Gland for Thermocouple



Thermocouple Entry Gland TCEG



Y Piece Extension



Tee Piece Extension



Clamps, Tri-clover Blanks, Gaskets



Type T thermocouple Wire



Steam Test Kits

Datasheet: 043

Thermal Detection Limited Unit 6, Primrose Hill Industrial Estate Orde Wingate Way Stockton-on-Tees, England, TS19 0GA



Fax: +44 (0) 1642 618307 Tel: +44 (0) 1642 602878 Email: sales@thermal-detection.com

Web: www.thermal-detection.com

Web: hermal Detection Ltd 2016